


Intended use . Kinetic system for Aspartate Aminotransferase (AST) or glutamic oxalacetic transaminase (GOT) determination.

Test principle . AST catalyzes specifically the transfer of the amino group of aspartic acid to a ketoglutarate yielding glutamate and oxalacetate. The oxalacetate is reduced to malate by the enzyme malate dehydrogenase (MDH), and coenzyme NADH is oxidized to NAD. The absorbance reduction at 340 nm, as a consequence of NADH oxidation, is determined photometrically and is direct proportional to the AST activity in the sample.

AST determination involves the following reactions:



Summary . The first U.V. kinetic methodology for measuring AST activity was described by Karmem¹ and improved later by Henry et al². This measuring procedure was widely accepted for being fast, precise and accurate.

The Reference Procedure proposed by International Federation of Clinical Chemistry and Laboratory Medicine (IFCC)³ recommends the use of pyridoxal phosphate to assure the total activation of all the transferase in the serum, avoiding false low results in samples with coenzyme deficiency. The Labtest system is improved following these recommendations.

In order to achieve IFCC Procedure traceable results, is needed the use of the two-reagent method with activation by the pyridoxal phosphate (Reagent 3). When the monoreagent is applied, the pyridoxal phosphate activation is not employed. Thus, the obtained results are not traceable to the IFCC Procedure.

The substances for this reaction are distributed properly in three reagents in order to get more stability in the original liquid form and keep the optimum conditions of operation conditions.

AST/GOT Labtest uses the kinetic system and is easily applied to most automatic and semi-automatic equipments which are able to measure the absorbance at 340 nm.

Methodology . U.V. Kinetic - IFCC

Reagents

1. [R1] - Reagent 1 - Store at 2 - 8 °C. Reagent label bears expiration date.

This reagent contains Tris buffer (105 mmol/L), L-aspartate (330 mmol/L), MDH (≥ 825 U/L), LDH (≥ 1200 U/L) and sodium azide (0.095%).

2. [R2] - Reagent 2 - Store at 2 - 8 °C. Reagent label bears expiration date.

This reagent contains Tris buffer (20 mmol/L), NADH (1320 μ mol/L), α -ketoglutarate (66 mmol/L) and sodium azide (0.095%).

3. [R3] - Reagent 3 - Store at 2 - 8 °C. Reagent label bears expiration date.

This reagent contains tris buffer (20 mmol/L), pyridoxal phosphate (11.1 mmol/L) and sodium azide (0.095%).

Precautions and warnings

For in vitro diagnostic use.

Disposal of all waste material should be in accordance with local guidelines.

The usual security cares should be applied on the reagent handling.

In automatic equipments, the reagents may be contaminated with other reagents or the air, depending on the equipment's characteristic and the work conditions. These can result in stability reduction and calibration modifications.

As it occurs in all enzymatic activity measurement, the incubation time and temperature is important for the quality of the results.

The reagents contain sodium azide as preservative. Avoid ingestion. In case of eyes contact, immediately flush eyes with plenty of water and get medical assistance.

Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide accumulation.

Storage and stability . Unopened reagents, when stored at indicated temperature, are stable up to expiration date shown on the label.

Deterioration . Microbial or chemical contamination may decrease reagents stability.

AST Work Reagent and Reagent 1, Reagent 2 and Reagent 3 mixture are not suitable for use if it has an absorbance lower than 1.0 at 340 nm when measured versus water as reference or in case of contaminations signs or if it develops turbidity.

Specimen collection and preparation

Use serum or plasma (EDTA, Heparin). AST is reportedly stable in serum or plasma (EDTA) for about 4 days at 2 - 8 °C and 2 weeks at -10 °C.

No known test method can offer complete assurance that human blood samples will not transmit infectious diseases. Therefore, all blood derivatives should be considered potentially infectious.

Interference

Chronic alcoholism increases AST activity in 18 to 100%. Anabolic steroids, chloramphenicol, chlorothiazide, prolonged use of aspirin, and others may increase the AST activity.

Bilirubin up to 19 mg/dL, hemoglobin up to 45 mg/dL, triglycerides up to 650 mg/dL and pyruvate up to 0.2 mmol/L do not interfere significantly. Values of hemoglobin greater than 45 mg/dL provide false increased results.

Hyperlipemic samples and jaundice remarkably increase the absorbance at 340 nm. When the enzyme activity of these kinds of samples is very elevated, the substrate consumption is very fast, without a significant decrease of the absorbance. In these samples, when a difference of absorbance per minute is very small, the determination should be repeated using a sample diluted with 0.85% NaCl.

Materials required not provided

1. Photometer capable of keeping the cuvette temperature at 37 °C, and measuring absorbance at 340 nm.
2. Pipetts to measure reagents and samples.
3. Timer.

Manual procedure

See **Calculation, Calibration, Measurement/ Reportable Range and Notes.**

Enzymatic Activity Determination using the pyridoxal phosphate In order to achieve traceable results to IFCC Procedure, is needed the use of the two-reagent method, to occur the enzyme total activation by the pyridoxal phosphate.

Preparing the reagent . Add 0.300 mL of the Reagent 3 to a bottle of Reagent 1 (24 mL) and mix gently. Stability: 21 days at 2 - 8 °C and 24 hours at 15 - 25 °C when no chemical or microbial contamination occurs. Optionally, a lower volume of the mixture (Reagent 1 + Reagent 3) may be prepared by using one part of the Reagent 3 to 80 parts of Reagent 1.

Procedure

1. In a test tube labeled "Test" or "Calibrator", add 0.800 mL of the mixture Reagent 1 + Reagent 3.
2. Add 0.100 mL of the sample or enzymes calibrator, homogenize and incubate in a water-bath at 37 ± 0.2 °C. Wait five minutes. After this incubation it is possible wait until 30 minutes to start the kinetic determination with the addition of the Reagent 2.
3. Perform a water blank measurement at 340 nm.

4. Add 0.200 mL of the Reagent 2, homogenize and transfer immediately to a cuvette at 37 ± 0.2 °C. Wait one minute.

5. Measure the initial absorbance (A_1), and start simultaneously the timer. Measure the absorbance again after 2 minutes (A_2). In order to verify the reaction linearity, it is recommended to measure in 1 minute as well, and check if the difference of absorbance in each minute is constant.

Enzymatic Activity Determination without the pyridoxal phosphate.

Preparing the working reagent . Transfer all the contents of one Reagent 2 bottle to one Reagent 1 bottle and mix gently.

The Working Reagent is stable 10 days at 2 - 8 °C and 24 hours at 15 - 25 °C, when no chemical or microbial contamination occurs. Optionally, a lower volume of the Working Reagent may be prepared by using the volume proportion 4:1 of the Reagent 1 and Reagent 2, respectively.

Procedure

1. In a test tube labeled "Test" or "Calibrator", add 1.0 mL of the Work Reagent.
2. Perform a water blank measurement at 340 nm.
3. Add 0.100 mL of the sample or enzymes calibrator, homogenize and transfer immediately to a cuvette at 37 ± 0.2 °C. Wait one minute.
4. Measure the initial absorbance (A_1), and start simultaneously the timer. Measure the absorbance again after 2 minutes (A_2).

In order to verify the reaction linearity, it is recommended to measure in 1 minute as well, and check if the difference of absorbance in each minute is constant.

Initial absorbance (A_1) equal to or lower than 0.8 indicates that the sample has elevated AST activity. In this case, dilute the sample and measure again (see Measurement/Reportable Range).

Calibration

Manual calibrations

Use Labtest Calibra series. AST activity is traceable to the IFCC reference method³.

Calibration frequency

Two or three point calibration after reagent lot change;
Two or three point calibration when the internal quality control indicates.

Automatic systems

Reagent Blank: reagent water or 0.85% NaCl.
Use Labtest Calibra series. AST activity is traceable to the IFCC reference method³.

Calibration frequency

Two or three point calibration after reagent lot change;
Two or three point calibration when the internal quality control indicates.

Quality control . For quality control use Qualitrol Level 1 and Qualitrol Level 2 or other suitable control material. The limits and control interval must be adapted to the laboratory requirements. Each laboratory should establish corrective actions to be taken if values fall outside the control limits.

Calculations . It is a usual procedure calculates the enzymatic activity results using a theoretical factor achieved in reaction optimum conditions, described below:

wavelength: 340 nm.
cuvette at 37 ± 0.2 °C, 10 mm light path.
pass band ≤ 2 nm.
stray light ≤ 0.1%.

If one of the correlated parameters is modified, it is recommended to apply an enzymes calibrator indicated by the reagent manufacturer. Labtest Diagnostica recommends Calibra series to perform the AST/GOT system calibration.

$$\Delta A/\text{minute (test or Calibrator)} = \frac{A_1 - A_2}{2}$$

$$\text{Factor} = \frac{\text{Calibrator activity}}{\Delta A/\text{minute Calibrator}}$$

$$\text{AST activity (U/L)} = \Delta A/\text{minute (test)} \times \text{Factor}$$

If all the correlated parameters are fulfilled, the theoretical factor (1746) can be applied.

Measurement/reportable range

From 3.5 up to 400 U/L.

If AST activity exceeds 400 U/L, the sample must be diluted with 0.85% NaCl. Multiply the result by the appropriate dilution factor.

Dilute the sample so that the obtained value is around 50 and 200 U/L.

Expected values⁷ . These values should be used as guidance only. Each laboratory should evaluate the transferability of the expected values to its own patient population and, if necessary, estimate its own reference interval.

Age	Men (U/L)	Women (U/L)
01-07 days	26-98	20-93
08-30 days	16-67	20-69
01-06 months	16-62	16-61
07-12 months	16-52	16-60
01-03 years	16-57	16-57
04-06 years	10-47	10-47
07-15 years	10-41	05-36
Adults	11-39	10-37

Conversion: Conventional Unit (U/L) x 16.7 = Unit IS (nkat/L).

Performance characteristics

Recovery studies . In two samples with aspartate aminotransferase activities of 162 and 261 U/L were added different quantities of the enzyme. Subsequent analyses provided recoveries ranging from 99.5 to 100%. The mean proportional systematic error at 77 U/L decision level was 0.2 U/L and at 285 U/L decision level was 0.7 U/L.

Method Comparison . The proposed method was compared with IFCC³ reference method, obtaining the follow results:

Tests performed using the pyridoxal phosphate

	Comparison Method	Labtest Method
N	40	40
Range (U/L)	17.4-296.6	18.8-287.4
Mean (U/L)	116.7	115.5
Regression analysis	Labtest Method = 0.995 x Comparison Method - 0.547	
Correlation coefficient	0.999	

The total systematic error (bias) is 1.2% in a decision level of 77 U/L and 0.7% in a decision level of 312 U/L.

Tests performed without the pyridoxal phosphate

	Comparison Method	Labtest Method
N	40	40
Range (U/L)	17.4-291.1	22.0-278.2
Mean (U/L)	113.5	109.5
Regression analysis	Labtest Method = 0.954 x Comparison Method + 1.147	
Correlation coefficient	0.997	

The total systematic error (bias) is 3.1% in a decision level of 77 U/L and 4.2% in a decision level of 285 U/L.

Imprecision

Tests performed using the pyridoxal phosphate

Imprecision - within run

	N	Mean (U/L)	SD (U/L)	(%) CV
Sample 1	20	77	1.71	2.2
Sample 2	20	312	2.69	0.9

Imprecision - run-to-run

	N	Mean (U/L)	SD (U/L)	(%) CV
Sample 1	20	77	1.58	2.1
Sample 2	20	312	6.85	2.2

Tests performed without the pyridoxal phosphate

Imprecision - within run

	N	Mean (U/L)	SD (U/L)	(%) CV
Sample 1	20	77	1.56	2.0
Sample 2	20	285	2.25	0.8

Imprecision - run-to-run

	N	Mean (U/L)	SD (U/L)	(%) CV
Sample 1	20	77	3.05	4.0
Sample 2	20	285	7.83	2.8

Analytical sensitivity . Detection limit: 1.75 U/L.
The detection limit represents the lowest measurable AST activity that can be distinguished from zero.

Effects of matrix dilution . Two samples with values 372 and 425 U/L were used for evaluating the system response on dilution with 0.85% NaCl. Using dilution factors on the range of 2 to 16, mean recovery was found to be 100.1 %.

Notes

- 1. The material cleaning and drying are fundamental factors to the reagent stability and to obtain correct results.
- 2. The water in the laboratory to prepare reagents and use in the measurements, must have resistivity ≥1 megaohm, or conductivity ≤1 microsiems and silicates concentration must be <0.1mg/L (Type II reagent water). The water for washing must be Type III, having resistivity ≥0.1 megaohms or conductivity ≤10 microsiemens. For the final washing, use Type II reagent water.
- 3. It is suggested to consult: <http://www.fxol.org> in order to review physiopathological source and drugs interference in results and methodology.

References

1. Karmen A. J Clin Invest 1955;34:131.

2. Henry RJ, Chiamori N, Golub O, Berkman S. Amer J Clin Path 1960; 34:381.

3. IFCC Reference Procedure for the Measurement of Catalytic Concentration of Aspartate Aminotransferase. Clin Chem Lab Med 2002; 40 (7):725-33.

4. Sociedad Española de Bioquímica Clínica y Patología Molecular, Base de Datos de Variación Biológica. Disponível em: <<http://www.seqc.es/article/articleview/330/1/170>> (access in 2006/04).

5. <http://www.westgard.com/biodatabase1.htm> (access in 2008/03).

6. Westgard JO, Barry PL, Hunt MR, Groth T. Clin Chem 1981; 27:493-501.

7. Soldin SJ, Brugnara C, Wong EC. Pediatric Reference Intervals, 5.ed. Washington: AACC Press, 2005. p. 39-40.

8. Labtest: data on file.

9. Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, 2.ed. Philadelphia: W.B. Saunders Company, 1994. p. 788-96.

Presentation

Product	Reference	Contents	
AST/GOT Liquiform	109-4/30	R 1	4 X 24 mL
		R 2	4 X 6 mL
		R 3	1 X 1,5 mL
	109-2/100	R 1	2 X 80 mL
		R 2	2 X 20 mL
		R 3	1 X 2,2 mL

Application procedures using AST Kinetic System are available for various automated instruments.

The number of tests in automated instruments depends on the programmed parameters.

Consumer information

[Warranty conditions]

Labtest Diagnóstica warrants the performance of this product under the specifications until the expiration date shown in the label since the application procedures and storage conditions, indicated on the label and in this insert, have been followed correctly.



Labtest Diagnóstica S.A.

CNPJ: 16.516.296 / 0001 - 38
Av. Paulo Ferreira da Costa, 600 - Vista Alegre - CEP 33400-000
Lagoa Santa - Minas Gerais Brasil - www.labtest.com.br

Consumer Service | e-mail: sac@labtest.com.br

Revision: May, 2013
Ref.: 180613

Copyright by Labtest Diagnóstica S.A.
Reproduction under previous authorization

Símbolos utilizados com produtos diagnósticos in vitro

Símbolos usados con productos diagnósticos in vitro

Symbols used with ivd devices

	Conteúdo suficiente para < n > testes Contenido suficiente para < n > tests Contains sufficient for < n > tests		Risco biológico Riesgo biológico Biological risk
	Data limite de utilização (aaaa-mm-dd ou mm/aaaa) Estable hasta (aaaa-mm-dd o mm/aaaa) Use by (yyyy-mm-dd or mm/yyyy)		Marca CE Marcado CE CE Mark
	Material Calibrador Material Calibrador Calibrator Material		Tóxico Tóxico Poison
	Material Calibrador Material Calibrador Calibrator Material		Reagente Reactivo Reagent
	Limite de temperatura (conservar a) Temperatura limite (conservar a) Temperature limitation (store at)		Fabricado por Elaborado por Manufactured by
	Representante Autorizado na Comunidade Europeia Representante autorizado en la Comunidad Europea Authorized Representative in the European Community		Número do lote Denominación de lote Batch code
	Consultar instruções de uso Consultar instrucciones de uso Consult instructions for use		Controle Control Control
	Número do catálogo Número de catálogo Catalog Number		Controle negativo Control negativo Negative control
	Adições ou alterações significativas Cambios o suplementos significativos Significant additions or changes		Controle positivo Control positivo Positive control
	Produto diagnóstico in vitro Dispositivo de diagnóstico in vitro In vitro diagnostic device		Controle Control Control
	Liofilizado Liofilizado Lyophilized		Corrosivo Corrosivo Corrosive
	Período após abertura Período post-abertura Period after-opening	Ref.: 201112	