COVID-19 Ag Rapid Test (Swab Naso/Oro)

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Ref.: 743

ANVISA 10009010374

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Referências

- World Health Organization. (2020). Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays: interim guidance, 11 September 2020. World Health Organization. <https://apps.who.int/iris/handle/10665/334253>. License: CC BY-NC-SA 3.0 IGO.
- Westgard JO, Barry PL, Hunt MR, Groth T. Clin Chem 1981;27:493-501.
- Labtest: Dados de Arquivo
- Cheng AC, Williamson DA. An outbreak of COVID-19 caused by a new coronavirus: what we know so far. Med J Aust 2020; 9 March.
- Xiang F, Wang X, He X et al. Antibody Detection and Dynamic Characteristics in Patients with COVID-19. Clin Infect Dis. 2020 Apr 19;ciaa461.
- Oliveira BA, Oliveira LC, Sabino EC, Okay TS. SARS-CoV-2 and the COVID-19 disease: a mini review on diagnostic methods. Rev. Inst. Med. trop. S. Paulo vol.62 São Paulo 2020 Epub June 29, 2020.
- Diao B, Wen K, Chen J et al. Diagnosis of Acute Respiratory Syndrome Coronavirus 2 Infection by Detection of Nucleocapsid Protein running title: Diagnosis of COVID-19 by N antigen detection. MedRxiv preprint doi: <https://doi.org/10.1101/2020.03.07.20032524>. March 13, 2020.
- Mathuria JP, YadavR, Rajkumarc. Laboratory diagnosis of SARS-CoV-2

- A review of current methods. J Infect Public Health. 2020 Jul; 13(7): 901–905.

9. Centers for Disease Control and Prevention. Overview of Testing for SARS-CoV-2 (COVID-19). <<https://www.cdc.gov/coronavirus/2019-ncov/hcp/testing-overview.html>>. Acessado em 02 de setembro de 2020.

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Referencias

1. World Health Organization. (2020). Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays: interim guidance, 11 September 2020. World Health Organization. <https://apps.who.int/iris/handle/10665/334253>. License: CC BY-NC-SA 3.0 IGO.
2. Westgard JO, Barry PL, Hunt MR, Groth T. Clin Chem 1981;27:493-501.
3. Labtest: Datos de Archivo
4. Cheng AC, Williamson DA. An outbreak of COVID-19 caused by a new coronavirus: what we know so far. Med J Aust 2020; 9 March.
5. Xiang F, Wang X, He X et al. Antibody Detection and Dynamic Characteristics in Patients with COVID-19. Clin Infect Dis. 2020 Apr 19;ciaa461.
6. Oliveira BA, Oliveira LC, Sabino EC, Okay TS. SARS-CoV-2 and the COVID-19 disease: a mini review on diagnostic methods. Rev. Inst. Med. trop. S. Paulo vol.62 São Paulo 2020 Epub June 29, 2020.
7. Diao B, Wen K, Chen J et al. Diagnosis of Acute Respiratory Syndrome Coronavirus 2 Infection by Detection of Nucleocapsid Protein running title: Diagnosis of COVID-19 by N antigen detection. MedRxiv preprint doi: <https://doi.org/10.1101/2020.03.07.20032524>. March 13, 2020.

8. Mathuria JP, YadavR, RajkumarC. Laboratory diagnosis of SARS-CoV-2 - A review of current methods. J Infect Public Health. 2020 Jul; 13(7): 901–905.

9. Centers for Disease Control and Prevention. Overview of Testing for SARS-CoV-2 (COVID-19). <<https://www.cdc.gov/coronavirus/2019-ncov/hcp/testing-overview.html>>. Acessado em 02 de setembro de 2020.

10. World Health Organization. Advice on the use of point-of-care immunodiagnostic tests for COVID-19. 2020.

11. Xu J, Zhao S, Teng T, et al. Systematic Comparison of Two Animal-to-Human Transmitted Human Coronaviruses: SARS-CoV-2 and SARS-CoV. Viruses 2020, 12, 244.

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Insert

Ref.: 743

Purpose . The product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743 is a rapid chromatographic immunoassay for the qualitative detection of SARS-CoV-2 antigens in nasopharyngeal or oropharyngeal swab samples as an aid in the diagnosis of new coronavirus infection.

Professional use only.

[For in vitro diagnostic use only].

Principle . The product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743 is a qualitative immunoassay for the detection of SARS-CoV-2 antigens in nasopharyngeal or human oropharyngeal swab samples. The product consists of a membrane containing anti-SARS-CoV-2 antibodies in the test line (T). Specimen from nasopharyngeal or oropharyngeal swabs react with anti-SARS-CoV-2 antibodies bound to colored particles that migrate by capillary action and form a colored line in the test (T) region if the sample contains virus antigens, indicating a positive result. If the sample does not contain SARS-CoV-2 antigens, no colored lines will appear in the test line region, indicating a negative result. As a procedural control, a colored line must always appear in the control line region (C).

System characteristics . The product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743 is an immunochromatographic method that allows the detection of SARS-CoV-2 antigens in nasopharyngeal and oropharyngeal swab samples.

In contrast to molecular and serological methodologies, the detection of antigens has the advantages of being a simple, quick and easy-to-interpret procedure, in addition to preceding the appearance of antibodies in infected people. Therefore, viral antigen detection can play a rapid screening effect and achieve the purpose of an early diagnosis of COVID-19.

Due to the advance of SARS-CoV-2 infections, the WHO started to recommend rapid tests for detection of the viral antigen, as they can play an important role in patient guidance and management, in public health decision-making and in the surveillance of COVID-19.¹

Methodology

Immunochemistry.

Reagents

1. Buffer – Store at 4 – 30 °C. Do not freeze it.

It contains Tris buffer, surfactant, stabilizer and denaturing agent. The buffer must only be opened at the time of use.

2. Test device – Store at 4 – 30°C. Do not freeze them.

They contain anti-SARS-CoV-2 antibody conjugated to colored microparticles as capture reagent and anti-SARS-CoV-2 antibody as detection reagent.

Once opened, the reagent plate is stable for up to one hour at room temperature.

Unopened reagents, when stored under the indicated conditions, are stable until the expiration date printed on the label. While handling, reagents are subject to chemical and microbial contamination that can cause reduced stability.

Storage and stability . The product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743 must be stored in hermetically sealed package at a temperature between 4 - 30 °C.

The product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743 is stable until the expiry date printed on the package. Do not use it after expiry date.

The product must remain in the sealed package until use. Once opened, the reagent plate (Ref. 743.2) remains stable for up to one hour at room temperature (15 - 30 °C). Discard the plate if the package is damaged.

The product must not be frozen or stored at a temperature above 30 °C. Humidity and temperature can adversely affect results.

Precautions and warnings

These Instructions for Use must be read completely before performing the test. The sample volume informed in the "Procedure" section must be strictly followed. Failure to follow the instructions for use may lead to incorrect results.

Since no known test can ensure that nasopharyngeal or oropharyngeal swab samples do not transmit infections, they should be considered potentially infectious. Thus, when handling them, biosafety rules must be followed.

To dispose reagents and biological material, we suggest applying local, state or federal protection standards. Used test must be disposed of in accordance with local, state, or federal environmental protection regulations.

The usual safety precautions should be used when handling reagents and samples. Wear protective clothing such as lab coats, disposable gloves and eye protection when samples are analyzed.

Required material and not provided

Chronometer.

Sample

Nasopharyngeal and Oropharyngeal Swab


Samples must be tested immediately after collection. If samples are not processed immediately, swabs should be placed in a dry, sterile, and tightly closed plastic tube for storage. Do not use viral transport media for swab storage. The swab sample is stable for up to 2 hours at temperatures between 18 - 28 °C or between 2 - 8 °C.

A Standard Operating Procedure (SOP) should be created to establish adequate procedures for sample collection, preparation and storage. We emphasize that the errors due to the sample can be much larger than the errors that occurred during the analytical procedure.

Sample collection


Nasopharyngeal swab

1. With the patient's head slightly tilted (45°), insert a sterile swab into the patient's nostril, reaching the surface of the posterior nasopharynx.
2. Rub the surface of the posterior nasopharynx for a few seconds to absorb nasal fluids.
3. Carefully remove the swab from the nasal cavity.


Oropharyngeal swab

1. With the head slightly tilted, the patient should open the mouth to reveal the pharyngeal tonsils.
2. Insert a sterile swab into the patient's oral cavity and rub it at least 3 times on the base of the tongue, each pharyngeal tonsil, and the posterior pharyngeal wall.
3. Carefully remove the swab from the oral cavity.

Sample preparation

1. Open the buffer tube and transfer the entire volume to the dropper bottle.
2. Insert the swab sample into the buffer and rotate it 10 times pressing it against the inner wall of the tube to release the antigens.
3. Remove the swab while squeezing the sides of the tube to extract all liquid from the swab.
4. Close the sample preparation tube.


Test procedure

The samples and the product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743 should be at room temperature (15 - 30 °C) at the time of running the test.

1. Remove the reagent plate from the sealed aluminum package and use it for a maximum period of one hour. The best results will be obtained if the test is performed immediately after opening the package.
2. Invert the dropper bottle containing the sample over the sample well and add 5 drops. Then, start the timer.

3. Read the result in 10 - 15 minutes. Do not interpret the result after 20 minutes.

Interpretation of results

See the illustration above for reference.

POSITIVE: Two distinct colored lines appear, one colored line in the control region (C) and one in the test region (T). The positive result can be confirmed as soon as the T and C lines appear indicating the detection of SARS-CoV-2 antigens in the sample.

The intensity of the color in the test line region (T) will vary according to the amount of SARS-CoV-2 antigen present in the sample. Therefore, any color tone in the test region (T) should be considered positive.

NEGATIVE: A colored line appears in the control region (C) and no visible colored lines appear in the test line region (T).

INVALID: There is no appearance of the control line. Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure to appear. Review the procedure and repeat the test with a new reagent plate.

Any sample that shows a positive result in a rapid test, due to this being a screening test, should be analyzed by complementary methodology, according to the flowcharts described by the World Health Organization.¹⁰

This product is a screening test, so a negative result does not rule out the possibility of SARS-CoV-2 infection. For the interpretation of the test result, other clinical information must be taken into account.

Interferences

No interference was noted for samples containing whole blood (1%) and mucin (3%).

No interference was observed for samples containing the drugs: nasal spray (20%), nasal corticosteroids (20%), Olopatadine Hydrochloride (20%), mouthwash (5 mg/mL), Tamiflu (5 mg/mL), Oseltamivir (5 mg/mL), Acetaminophen (10 mg/mL), Ibuprofen (10 mg/mL), Povidone-iodine (1%), Acetylsalicylic acid (20 mg/mL) and Cefadroxil (5 mg/mL).

Test limitations

1. Test performance depends on the quality of the sample. Failure to follow the procedure may lead to incorrect results.

2. The product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743 is intended only as an aid to in vitro diagnostics. This test should be used for detection of SARS-CoV-2 antigens in nasopharyngeal or oropharyngeal swab samples as an aid in diagnosing patients with suspected SARS-CoV-2 infection in conjunction with clinical findings and results of other laboratory tests. Neither the quantitative value nor the increase in concentration of SARS-CoV-2 antigens can be determined by this test.

3. The product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743 indicates only the presence of SARS-CoV-2 antigens in the sample and should not be used as the sole criterion for the diagnosis of infections caused by SARS-CoV-2.

4. Negative results do not exclude SARS-CoV-2 infection, particularly in those individuals who have had contact with the virus. If the test result is negative or non-reactive and clinical symptoms persist, it is recommended to collect a new sample from the patient for retesting. Follow-up tests with molecular diagnosis should be considered to rule out infection in these individuals.

5. The test may give negative results if the titre of the SARS-CoV-2 antigens in the sample is below the detection limit of the test (10 ng/mL of recombinant SARS-CoV-2 protein or viral titre of $6.12 \times 10^{2.2}$ TCID₅₀/mL).

6. Infections with non-SARS-CoV-2 variants of coronaviruses may cause false positive results.¹¹

Internal quality control². The laboratory must maintain an internal quality control program that clearly defines the applicable regulations, objectives, procedures, criteria for quality specifications and tolerance limits, corrective actions and record of activities.

The formation of a colored line in the control region is indicative of adequate performance of the procedure. To ensure that the system has not been adversely affected and maintains the established performance

levels, we suggest associating a quality control system using, daily, known samples, one positive and one negative.

Performance characteristics³

Comparative Studies

The product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743 was evaluated with 734 samples being 249 positive and 485 negative for SARS-CoV-2.

COVID-19 Ag Rapid Test (Swab Naso/Oro) Ref. 743	RT-PCR		
	Positive	Negative	Total
	Positive	228	9
Negative	21	476	497
Total	249	485	734

Sensitivity: 91.6% (IC 95%: 87.4% to 94.7%)

Specificity: 98.1% (IC 95%: 96.5% to 99.2%)

Accuracy 95.9% (IC 95%: 94.2% to 97.2%)

Repeatability - intra-assay inaccuracy . Intraassay imprecision was verified by evaluating ten replicates of two different samples, one negative and one positive for the SARS-CoV-2 antigen. The trials were carried out on three batches of the product in accordance with its instructions for use. Negative and positive results found demonstrate perfect agreement with the expected results.

Reproducibility - Inter-assay inaccuracy . The total inaccuracy was verified in three tests carried out on consecutive days, totaling 30 replicates, using three different batches of the product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743. Ten trials per day were carried out, evaluating two different samples, one negative and one positive for the SARS-CoV-2 antigen. The trials were carried out in accordance with the instructions for use of the product. Negative and positive results found demonstrate perfect agreement with the expected results.

Pro-Zone Effect . The product COVID-19 Ag Rapid Test (Swab Naso/Oro) - Ref. 743 has no prozone effect in samples with SARS-CoV-2 titer up to $1 \times 10^{6.2}$ TCID₅₀/mL.

Detection limit . The detection limit of the product COVID-19 Ag Rapid Test (Swab Naso/Oro) Ref. 743 is 10 ng/ml for recombinant SARS-CoV-2 protein or viral titer of $6.12 \times 10^{2.2}$ TCID₅₀/mL.

Cross Reactivity . No cross-reaction was noted when samples enriched with *Chlamydia pneumoniae*, *Streptococcus pyrogenes*, *Mycoplasma pneumoniae*, *Haemophilus influenzae*, *Pseudomonas aeruginosa*, *Mycobacterium tuberculosis*, *Legionella pneumophila*, *Staphylococcus epidermis*, *Bordetella pertussis*, *Streptococcus pneumoniae*, *Staphylococcus salivarius*, *Candida albicans*, Coronavirus - FCV(3A2), Coronavirus - FIP(2A4), Influenza A H3N2, Influenza B, Coxakie A Virus, Coxakie A4 Virus, Coxakie B Virus, Coxakie B Virus, Adenovirus type 1, Adenovirus type 7, Enterovirus 71, Rhinovirus RV21, Respiratory syncytial virus A, Echovirus 6, Echovirus 9, Human Cytomegalovirus, Herpesvirus 1, Herpesvirus 2, Measles Virus, Mumps

Virus and Poliovirus – sabin (3A4) were tested with the COVID-19 Ag Rapid Test product (Swab Naso/Oro) - Ref. 743.

Clinical significance . The new coronavirus, formally named severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), is the causative agent of the 2019 coronavirus disease (COVID-19) that emerged in Wuhan City, Hubei Province, in Central China, December 2019.⁴ Molecular diagnosis, RT-PCR, is the gold standard for confirming SARS-CoV-2 infection in samples from the upper respiratory tract (nasopharyngeal and oropharyngeal secretions). As a complementary approach, IgG and IgM antibodies against SARS-CoV-2 can be detected in the intermediate and later stages of the disease (fourth day after the onset of symptoms).^{5,6} The biggest problem faced in the material-based diagnostic model Viral genomic is associated with high operational cost and the need for highly specialized labor in carrying out the trials. In the case of immunological tests, the limitations are their use in the initial phase of the infection because at that time the immune response is still in formation.⁷ In contrast to the methodologies presented, the detection of coronavirus antigen has the advantage of being fast and convenient, since the viral antigen is the specific marker of the virus and precedes the appearance of antibodies in infected people. Therefore, viral antigen detection may play a rapid screening effect and achieve the purpose of an early diagnosis of COVID-19.^{8,9} Rapid antigen tests are commonly used in the diagnosis of respiratory pathogens, including influenza virus and respiratory syncytial virus. The clinical performance of rapid diagnostic tests based on antigen detection depends on the circumstances in which they are used and they perform best when the person is tested in the early stages of SARS-CoV-2 infection when the viral load is generally higher. They can be used as screening tests in high-risk environments and quickly identify people infected with SARS-CoV-2, allowing for quick decision-making on infection prevention and control measures.⁹ It is important to emphasize that any sample that shows a positive result in a rapid test, due to this being a screening test, should be analyzed by complementary methodology, according to the flowcharts described by the World Health Organization.¹⁰

Presentation

Product	Referene	Content
COVID-19 Ag Rapid Test (Swab Naso/ Oro) Ref. 743	743-1	Buffer 1 x 0,35 mL
		Test Device 1 units
		swab 1 units
		Dropper bottle 1 units
	743-5	Buffer 5 x 0,35 mL
		Test Device 5 units
		swab 5 units
		Dropper bottle 5 units
	743-10	Buffer 10 x 0,35 mL
		Test Device 10 units
		swab 10 units
		Dropper bottle 10 units
	743-20	Buffer 20 x 0,35 mL
		Test Device 20 units
		swab 20 units
		Dropper bottle 20 units
	743-40	Buffer 40 x 0,35 mL
		Test Device 40 units
		swab 40 units
		Dropper bottle 40 units
743-100	Buffer 100 x 0,35 mL	
	Test Device 100 units	
	swab 100 units	
	Dropper bottle 100 units	

References

1. World Health Organization. (2020). Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays: interim guidance, 11 September 2020. World Health Organization. <https://apps.who.int/iris/handle/10665/334253>. License: CC BY-NC-SA 3.0 IGO.
2. Westgard JO, Barry PL, Hunt MR, Groth T. Clin Chem 1981;27:493-501.
3. Labtest: Dados de Arquivo
4. Cheng AC, Williamson DA. An outbreak of COVID-19 caused by a new coronavirus: what we know so far. Med J Aust 2020; 9 March.
5. Xiang F, Wang X, He X et al. Antibody Detection and Dynamic Characteristics in Patients with COVID-19. Clin Infect Dis. 2020 Apr 19;ciaa461.
6. Oliveira BA, Oliveira LC, Sabino EC, Okay TS. SARS-CoV-2 and the COVID-19 disease: a mini review on diagnostic methods. Rev. Inst. Med. trop. S. Paulo vol.62 São Paulo 2020 Epub June 29, 2020.
7. Diao B, Wen K, Chen J et al. Diagnosis of Acute Respiratory Syndrome Coronavirus 2 Infection by Detection of Nucleocapsid Protein running title: Diagnosis of COVID-19 by N antigen detection. MedRxiv preprint doi: <https://doi.org/10.1101/2020.03.07.20032524>. March 13, 2020.

8. Mathuria JP, YadavR, Rajkumarc. Laboratory diagnosis of SARS-CoV-2 - A review of current methods. J Infect Public Health. 2020 Jul; 13(7): 901–905.

9. Centers for Disease Control and Prevention. Overview of Testing for SARS-CoV-2 (COVID-19). <<https://www.cdc.gov/coronavirus/2019-ncov/hcp/testing-overview.html>>. Acessado em 02 de setembro de 2020.

10. World Health Organization. Advice on the use of point-of-care immunodiagnostic tests for COVID-19. 2020.

11. Xu J, Zhao S, Teng T, et al. Systematic Comparison of Two Animal-to-Human Transmitted Human Coronaviruses: SARS-CoV-2 and SARS-CoV. Viruses 2020, 12, 244.

Consumer information

[Guarantee Terms and Conditions]

Labtest Diagnóstica guarantees the product performance, within the specifications, until the expiration date indicated on the labels, provided that the care for use and storage indicated on the labels and in these instructions are correctly followed.

Labtest Diagnóstica S.A.

CNPJ: 16.516.296 / 0001 - 38

Av. Paulo Ferreira da Costa, 600 - Vista Alegre - CEP 33.240-152

Lagoa Santa - Minas Gerais Brasil - www.labtest.com.br

Customer Service | e-mail: customerservice@labtest.com.br

Edition: June, 2021

Revision: -

Ref.: 220621(00)

Copyright by Labtest Diagnóstica S.A.

Reproduction under previous authorization

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